PLS

PLS Hip
Stelo PLS

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History and Philosophy

Storia e filosofia

Since it was launched on the market in 1985, the finned stem designed by Prof. Spotorno has proven to be one of the most successful cementless prosthetic stems.

In 2003 Lima-Lto gave a new vision to the design creating the PLS stem, seeking the best proximal fixation without bone-implant contact in the distal portion of the stem to prevent stress-shielding.

The straight cementless stem features a rectangular section with sloping profiles which, in the proximal portion, presents fins longitudinal to the diaphysis that taper distally. Penetrating into the bony trabecular structure of the femur, these fins ensure an optimal primary stability.

The stem was also designed to minimise and prevent undesired rotation, varus and subsidence.

To this end, square-shaped sections, medial curvature and metaphyseal fins have been studied exhaustively in order to improve and optimise their design, creating a stem with excellent stability, adaptable to the various types of femoral geometry.

Moreover, its macro-roughened surface allows the bony integration of the stem, thus permitting an optimal secondary stability.

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Macro Rough surface finish
Finitura al corindone

The stem is manufactured in Ti6Al4V titanium alloy, according to strictly standardised production criteria and improved by the high-level technological know-how of Lima-Lto.

A preliminary phase of drop forging followed by a special thermal treatment make it possible to obtain the high mechanical performances required by the stem.

The subsequent mechanical machining by means of modern numerical control operation centres determines the final self-locking geometry that characterises the well-known bio-mechanical behaviour.

The external surface features a macro-roughened surface obtained through the projection of a controlled flow of corundum particles. Pressure, speed and granulometry of the casting make it possible to obtain a uniform and calibrated surface texture with appropriate rugosity values ($Ra\ 3\div 4\mu m$) for osteointegration.

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Biomechanical Characteristics

Caratteristiche biomeccaniche

The stem geometry and finishing guarantee an excellent primary stability as well as an uniform and lasting secondary osteointegration.

Three different anatomical configurations are available: "Standard", "Lateralising" and "Medialising". These configurations allow the surgeon to reconstruct the different joint anatomies by transferring the rotation centre from a "Standard" position (reference offset 39.5 to 55.3mm) to a "Lateralised" (offset +5mm) and to a "Medialised" (offset -5mm) position.

The possibility of increasing the offset without changing the whole length of the stem ensures optimal stability to the implant.

The wide medial curvature is adequately rounded for anatomical adaptation of the prosthesis to the calcar area. The distal end is further tapered in order to avoid any anomalous pressure cones.

An excellent biological fixation is reached thanks to the macro-roughened finishing of the surface.

The neck polished finishing reduces the probability of polyethylene particulate formation during accidental contact with the acetabular component. The neck is also tapered and, together with the 12/14 low coupling taper, allows the joint's range of motion to exceed 150° (for couplings with 36mm cephalic diameter)

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Applications

Cementless hip prostheses with biological fixation for the partial or total replacement of the hip where primary stability is guaranteed by means of the rectangular-section wedge-shaped longitudinal profile.

Ideal for the treatment of common hip pathologies.


Effective also in applications for cases of abnormal bony morphologies, whether primary (hip dysplasia) or secondary (trauma or osteotomy outcomes) morphologies, for all bony structures with adequate bony trophism.

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Sizes

The stems are available in 3 different configurations:

standard (132°), lateralising (offset +5mm) and medialising (offset -5mm).

For each of these configurations 11 progressively increasing sizes are available.

This permits to obtain a good functional result on all occasions, thanks to an optimal anatomical fitting of the prosthesis and to a balanced recovery of the joint's motions.

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
For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

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X Lima Deutschland GmbH