Cemented cups + Acetabular Plates
Coppe cementate + Gabbie

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
Standard Cup
Coppa standard

Protruded Cup
Coppa protrusa

The cemented acetabular cup are entirely produced in UHMWPE and are available with an outer diameter of 44, 46, 48, 50, 52, 54, 56, 58 mm. The inner diameters are appropriate to be coupled to femoral heads of 28 and 32 mm.

The acetabular cups are available in two different versions:

- ▶ standard;
- ▶ with an anti-dislocation shoulder, also for 36 mm diameter heads.

The standard version is completely hemispherical and is therefore indicated for deep acetabula.

The anti-dislocation version has a lateral portion which guarantees a greater coverage of the femoral head and prevents any accidental dislocation.

The identification diameters refer to the actual diameters of the cups.

To prepare the acetabular cavity it is therefore necessary to use reamers of larger diameters to ensure a uniform distribution of the cement mantle. It is possible to use the round reamers of the SPH series, whose identification diameter coincides with the actual diameter of the reamers themselves.

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The acetabular plates have been designed as supporting structures when cementing cups entirely made of UHMWPE.

Such reinforcing cages are made of pure Titanium with 50, 54, 58, 62 and 66 mm external diameters.

The possibility of using a plate in Titanium has been conceived to obviate the cases of intolerance to certain elements composing Steels (as Nickel, Chromium) and avoid any corrosion phenomenon.

The main indications for using the reinforcing cages are:

- ▶ revision of previous implants in presence of insufficient trophism or seriously altered bone structures;
- ▶ congenital dysplasias;
- ▶ acetabular fractures;
- ▶ in general, in any cases providing the use of cement for fixing the cup with reconstruction of the acetabular seat.

The reinforcing cages are anchoring to the pelvis by a caudal hook, to be inserted into the foramen obturatoris, two screws by the superolateral area and other screws fixed to the iliac wing with three cranial wings.

The wings are opportunely deformed during the operation to obtain the highest congruence between the implant and the anchoring surfaces.

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
The kind and length of screws better adapting to the considered case may be employed with the reinforcing cages. The use of Titanium screws is anyway recommended to prevent any possible corrosion phenomenon between the heads of the screws and the cage.

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


**Single Use Acetabular Cups
Similcotili monouso**

ref.		description	descrizione	qt.
9632.74.010	A	single use acetabular cup - dia. 50 mm	similcotile monouso - dia. 50 mm	1
9632.74.020	A	single use acetabular cup - dia. 54 mm	similcotile monouso - dia. 54 mm	1
9632.74.030	A	single use acetabular cup - dia. 58 mm	similcotile monouso - dia. 58 mm	1
9632.74.040	A	single use acetabular cup - dia. 62 mm	similcotile monouso - dia. 62 mm	1
9632.74.050	A	single use acetabular cup - dia. 66 mm	similcotile monouso - dia. 66 mm	1

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Operating Technique “Cemented Cups”

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(1)

1. Incision

There are many viable positions, but lateral decubitus is the one recommended for the MUELLER cups system.

The choice of an anterolateral access route places the incision – which stretches both distally and laterally to ensure adequate exposure of the articulation - laterally and proximally to the greater trochanter.

2. Acetabular Reaming

Expose the acetabulum with the surgical technique usually adopted to ensure an adequate view of the same for precise boring.

Remove former implants, if any.

Then remove any osteophytes and soft tissue surrounding the cotyloid rim to ensure an excellent view of bone anatomy and to check whether there are any cavity and/ or segmental defects.

If there are any bone defects, fix the bone graft before reaming the acetabulum.

Prepare the site of the cotyloid with the usual acetabulum reamer (Fig. 1).

Start reaming by choosing a reamer with an appropriate diameter. A reamer with a smaller diameter than the acetabulum's – which is defined

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cod. 9057.20.811

[2]

by means of transparencies – is to be preferred.

Follow the steps described to insert the spherical reamer into the reamer holder (Fig. 2):

- ▶ insert the reamer into the pegs present on the locking lever;
- ▶ push the reamer handle's locking track downwards and perform a quarter turn anti-clockwise with the reamer (Fig.3);
- ▶ then, release the locking track (Fig. 4)

[3]

Insert the reamer into the acetabulum maintaining, as far as possible, an abduction axis of approx. 45° and an anteversion axis of 15°. Move the reamer holder in a clockwise direction (Fig. 5).

[4]

For removal, lower the bayonet mount's locking ring nut, rotate the reamer clockwise and remove it from the holder.

Shift gradually to reamers with larger diameters to obtain an even spherical bone surface.

[5]

Note. Since the reamer is spherical, it can be rotated around any axis while in use to improve bone site preparation. Anyhow avoid antero-posterior or lateromedial traverse movements that could produce a

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non spherical acetabular conformation.

3. Testing and inserting the cup

Once the acetabular site's preparation phase has been completed, screw trial cup the size of the last reamer used onto the multipurpose handle or the positioner (Fig. 6). Use the wide openings available to ensure even contact along the acetabular rim.

The implanted cup will have the same nominal diameter as the last reamer used (Fig. 7).

Prepare surgical cement and position it in the acetabulum. Manually insert the cup and press it into the cement bed, taking care to protect the polyethylene articular surface (Fig. 8).

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cod. 9056.10.010

[9]

Use the cup impactor (Figs. 9-10), which has the same diameter as the head to be used, to apply pressure on the cup until cement polymerisation is complete. (Fig. 11)

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Tecnica operatoria
“Gabbie”

1. Acetabulum Reaming

Expose the acetabulum with the surgical technique usually adopted to ensure an adequate view of the same for precise boring.

Remove former implants, if any.

Then remove any osteophytes and soft tissue surrounding the cotyloid rim to ensure an excellent view of bone anatomy and to check whether there are any cavity and/ or segmental defects. In fact, correct positioning of the caudal hook requires clear exposure of the lower edge of the acetabulum.

If there are any bone defects, fix the bone graft before reaming the acetabulum.

Reamer with Cross Coupling
Fresa con attacco a croce

cod. 9057.20.942

[12]

Prepare the site of the cotyloid with the usual acetabulum reamer. (Fig. 12). Start reaming by choosing a reamer with an appropriate diameter. A reamer with a smaller diameter than the acetabulum's – which is defined by means of transparencies – is to be preferred.

Reamer Handle EZ
Manipolo per fresa EZ

cod. 9057.20.811

[13]

Follow the steps described to insert the spherical reamer into the reamer holder (Fig. 13):

- ▶ insert the reamer into the pegs present on the locking lever;

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[14]

- ▶ push the reamer holder handle's locking track downwards and perform a quarter turn anti-clockwise with the reamer (Fig. 14);
- ▶ then, release the locking track (Fig. 15).

[15]

Insert the reamer into the acetabulum maintaining, as far as possible, an abduction axis of approx. 45° and an anteversion axis of 15°. Move the reamer holder in a clockwise direction (Fig. 16).

For removal, lower the bayonet mount's locking ring nut, rotate the reamer clockwise and remove it from the holder.

Shift gradually to reamers with larger diameters to obtain an even spherical bone surface.

[16]

Note. Since the reamer is spherical, it can be rotated around any axis while in use to improve bone site preparation. Anyhow avoid anteroposterior or lateromedial traverse movements that could produce a non spherical acetabular conformation.

2. Testing and Inserting the Cage

Once the acetabular site's preparation phase has been completed, screw a trial cup the size of the last reamer used onto the multipurpose handle or

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[17]

Single use acetabular cup
Similcotile monouso

cod. 9632.74.010 ...

the positioner. Use the wide openings available to ensure even contact along the acetabular rim.

Then implant a cage the same nominal diameter as the last reamer used.

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[18]

At this point position the template cage (Fig. 17), having the same size as the cage to be implanted, inside the acetabular cavity.

Deform the wings and the caudal hook to suit acetabular morphology (Figs. 18-19).

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[22]

[23]

Remove the cage from the sterile package, choosing the same size as the last reamer used. Screw the multipurpose handle onto the cage's polar hole. Reproduce the template cup shape by folding the wings of the cage to be implanted with the handle for winglet modeling (Fig. 20) (pay attention not to deform the wings repeatedly to avoid them from strain hardening and, subsequently, breaking) (Figs. 21-22).

Enlarge the caudal hook to avoid over binding the system, which could cause strain breakage.

Apply pressure on the cage after filling any bone gaps with a bone graft (Fig. 23). Strike the placer axially by inserting the cage stably, ensuring that the caudal hook is correctly housed in the foramen obturatoris.

Check that the cup has adequately sunk into the acetabular site and test the implant's primary stability by applying gentle lever stress on the various planes with the placer rod.

Unscrew both placer and adaptor to ensure that there is contact between the cup and the bottom of the acetabulum through the bores present. If necessary, screw the rod on once again and press the cup further.

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cod. 9055.17.010 **[25]**

3. Bone Screw Fixation

Prepare the holes required to fix the cranial wings with the helix drills and the drill guide (Fig. 24).

To correctly insert the screws, position the track by points into the internal holes and perforate with the special points.

Use the number of screws that will guarantee correct mechanical stability. All holes can be used as long as there is an adequate bone stock to anchor the screws. To optimise cranial wing fixation, first link the proximal holes and then the distal ones.

Take each spongy screw with the special pincers and start screwing with the screwdrivers (Fig. 25).

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[26]

When the screw proceeds stably into the site, remove the pincers without locking it completely. Repeat the procedure for every screw you wish to insert.

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Perform final locking when all screws required have been inserted (Fig. 26).

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[27]

If the screws fitted into the cranial wings do not suffice to ensure correct primary stability, part of the holes inside the cage can be used to insert the spongiuous screws (Figs. 27-28).

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[28]

The internal holes have the primary function of allowing cement infiltration into the acetabular bone or into the bone graft.

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After fixing the cage and checking its stability (Fig. 29), cement the polyethylene cup by following the operating technique adopted for cemented cups.

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[29]

For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

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X Lima Deutschland GmbH
X Lima Austria GmbH