Fit

Anatomical Femoral Stem
Stelo femorale anatomico

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History and Philosophy

Storia e filosofia

The FIT cementless anatomical hip prosthesis has been designed based on the anatomical study of human femurs supported by experiences made in Germany and the results achieved using products already available on the market.


The prosthesis has been designed to obtain maximum stability while preserving femoral morphology and avoiding to force the implant component into the bone seat.

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Materials

Materiali

► The stem is made of Ti6Al4V titanium alloy.

► The meta-epiphyseal area is coated with porous titanium (approx. 400 μm), applied with vacuum spray plasma and a second thin layer of hydroxyapatite (60 μm).

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Biomechanical Characteristics

Caratteristiche biomeccaniche

► The prosthesis is placed on the meta-epiphysis depending to the femur curvature so as to ensure maximum torsion stability.


► The proximal articular load is distributed simultaneously both medially and laterally. The prosthesis is supported medially thanks to the anatomical curvature, while laterally it is supported by a marked flare (9°) at the base of the isthmus.

► An optimal centring of the stem is achieved. Thanks to the implant asymmetry, the distal part of the stem ensures a proper filling of the canal, at the same time preventing an anomalous load transfer to the bone structure and, therefore, undesired by-pass occurrences.

► The lowered and low tapered profile neck makes it possible to avoid accidental contact with the acetabular cup, and to increase the range of motion up to 144° with a head of 28mm. On the 12/14 taper, ceramic or metal modular heads of various lengths and diameters can be used.

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► The distal tract of the stem is highly polished in order to prevent any bony adhesion that would cause load by-pass to the proximal support.

► The neck is highly polished to minimize cup polyethylene abrasion in the event of abnormal contacts.

► The broaches provided with the specific instrument set are 5mm longer than the stems, so that the tapered stem tip is placed in line with the diaphysis and is free from any distal contact.

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
Stem
Stelo

Product Characteristics

- ▶ The FIT hip prosthesis system comprises a set of anatomical stems in progressively increasing sizes.
- ▶ The cervico-diaphyseal angle of 132° is the best for the most common cases of varus and valgus correction of the femoral neck.
- ▶ Thanks to the increase in neck size in relation to the stem size, it is possible to easily restore the femoral anatomic off-set.

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Indications

- ▶ Cementless anatomical prosthesis for total hip replacement
- ▶ Coxarthrosis in bone structures of adequate resistance

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Sizes

- ▶ 7 right sizes and 7 left sizes, increasing progressively
- ▶ Cervico-diaphyseal angle of 132°

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



Uncemented Porous-Coated Anatomic Total Hip Replacement. Results at Six Years in a Consecutive Serie.

Kim Yh - Kim Vem

(The Journal of Bone and Joint Surgery, 75-B: 6, 1993)

The Porous Coated Anatomic Total Hip Replacement: A Ten to Fourteen-Year Follow-up Study of a Cementless Total Hip Arthroplasty.

Haruo Kawamura - Michael J. Dunbar - Paraic Murray - Robert B. Bourne - Cecil H. Rorabeck

(The Journal of Bone and Joint Surgery, Am 2001 83: 1333-1338)

The Porous-Coated Anatomic Total Hip Prosthesis, Inserted without Cement. A Prospective Study with a Minimum Ten Years of Follow-up.

John S. Xenos - John J. Callaghan - R. David Heekin - William J. Hopkinson - Carlton G. Savory - Milan S. Moore

(The Journal of Bone and Joint Surgery, Am 1999 81: 74-82)

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P. Dalla Pria

S.I.O.T. 2002 - Lido di Venezia, Italy - 22/23 October 2002

Concept, Indications, Surgical Technique of the 'Fit' Anatomical Stem.

U. De Nicola

S.I.O.T. 2002 - Lido di Venezia, Italy - 22/23 October 2002

The 'Fit' Anatomical Stem. First Clinical and Radiological Results.

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