# SMR

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Cemented Revision Stem  
*Stelo da revisione cementato*

**Cemented Conical Stems**  
***Steli conici cementati***

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Uncemented Glenoids  
*Glenoidi non cementate*

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■ upon request / *su richiesta*

Cemented Glenoid  
*Glenoide cementata*


**Cemented Glenoids**  
***Glenoidi cementate***

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





Glenoids + Reverse

**RANDELLI REVERSE REVI-  
SION:**Glenoids + Reverse + Randelli Re-  
verse**ENDOPROTESI:**

Comune + Endoprotesi

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Comune + Endoprotesi + Glenoidi

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Comune + Glenoidi + Inversa

**REVISIONE SMR:**Comune + Revisione + Endoprotesi +  
Glenoidi + Inversa**INVERSA REVISIONE SMR:**

Glenoidi + Inversa

**INVERSA REVISIONE RAN-  
DELLI:**

Glenoidi + Inversa + Randelli Inversa

## “SMR Trauma” Surgical Technique

## Tecnica operatoria “SMR trauma”

The prosthetic treatment of acute traumas is generally limited to complex fractures/dislocations: split fractures with 3-4 fracture fragments, dislocation fractures of the humeral head, compression fractures with the destruction of more than 50% of the articular surface, etc.

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### 1. SURGICAL FIELD PREPARATION

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#### 1.1. Positioning

The patient is placed supine in the semi-sitting (beach chair) position (Fig. 1). The patient's arm must protrude out of the surgical table to allow a good mobility in all the planes. The patient's head must be fixed rigidly.

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#### 1.2. Access routes

Anterior, deltopectoral approach. The incision must be rather long, from the edge of the clavicle above the caracoid process, diagonally down to the humeral insertion of the deltoid (Fig. 1).

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The deltopectoral groove is localised and the cephalic vein isolated (moved medially or ligated and then dissected). After having incised the calvipectoral fascia, the conjoined tendon of the short head of the biceps and the coracobrachialis is moved medially.

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#### 1.3. Removal of fracture fragments

The presence of blood and lacerations of soft parts may cause some difficulties at this stage.

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(1)

It is advisable to localize distally the tendon of the long head of the biceps and then follow it proximally up to the rotator interval: the large tuberosity is localised laterally to the bicipital groove, while the small tuberosity is localised medially. These fragments should be carefully preserved with the inserted cuff tendons. The articular cavity is penetrated and the other fragments, as well as the humeral head, are removed. The arm may now be dislocated by external rotation and extension, allowing a full exposure of the diaphyseal stump.

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[2]

## **2. HUMERAL CANAL PREPARATION**

### **2.1a. Insertion of the trial stem with the beater-impactor (option 1)**

Screw the smaller diameter conical trial stem (B1) on the impactor (C1) (Fig. 2).

Then impact the stem into the canal. Next to the point where the stem fins are no more visible, insert the centimeter (T2) in the correspondent seat on the impactor handle and check the insertion depth (Fig. 3).

Then pass to the trial tapered stems with larger diameters until the correct fit inside the canal is obtained and the gauge reaches the required implant height.

This assessment is only indicative and does not determine the implant final height, which will be established in a second time with trial reductions.

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[3]

**[4]**

Making reference to a precise point of the humeral stump, bear in mind the measure indicated by the gauge. Once the correct diameter has been assessed, remove the beater-impactor leaving the stem in place. Insert and fix a **trial humeral body (G2)** - medium, long or short - with the correct retroversion, tightening the screw with the **allen wrench (M2)** (Fig. 4). Use the **body stopper (H2)** inserted in the cephalic cone, so as to avoid the humeral torsion (Fig. 5).

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[6]

[7]

### 2.1b. Insertion of the trial implant with the prosthesis introducer (option 2)

Place a trial medium humeral body (G2) on the smaller diameter conical trial stem (B1). Fix the two elements with the hexagonal wrench (Fig. 6).

Insert the trial implant into the prosthesis introducer (A2) and fix it by tightening the proximal screw (Fig. 7).

Before introducing the assembly in the canal, insert the alignment rod (P2) into the prosthesis introducer.

With the arm flexed at 90°, place the rod parallel to the forearm. In this position the implant will have a 30° retroversion (Fig. 8). If a smaller retroversion is required, the rod must be extra-rotated.

Impact the stem into the canal. Next to point where the stem fins are no more visible, insert the centimeter (T2) in the corresponding seat on the introducer handle (Fig. 9).

Then pass to the trial stems with larger diameters until the correct fit inside the canal is obtained and the gauge reaches the required implant height.

This assessment is only indicative and does not determine the implant final height, which will be established in a second time. Making reference to a precise point of the humeral stump, bear in mind the measure indicated by the gauge marks.

Once the trial stem with the correct diameter has been identified, remove the prosthesis introducer, leaving the stem and the trial body in place.

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
[10]

### 3. TRIAL ADAPTORS AND HUMERAL HEADS

Drive the **trial adaptor, STD (E2)** (neutral or eccentric) or **long** (neutral or eccentric) type, into the **trial head (C2)** by hand (Fig. 10). Insert the head in the cephalic cone of the trial body (Fig. 11).

Reduce the implant and check the contact with the glenoid.

The diameter of the head is determined by placing the **head gauge (D2)** on the resected bony head (Fig. 12).

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
[12]

[ 13 ]

In case an eccentric adaptor is used, mark the head position with an electric scalpel, using the arrow marked on the head as a reference.

This procedure allows the head to be placed in its final correct position (Fig. 13).

If there is no bone in the humerus near the head, bear in mind the last seat of the adaptor peg in the trial humeral body (Fig. 14). The peg of the final head must be reinserted in the same seat.

Carry out the trial reductions.

Use the appropriate **pliers (F2)** to remove the trial adaptors (Fig. 15).

[ 14 ]

#### **4. TRIAL REDUCTIONS AND MODIFICATION OF THE COMPONENTS**

Any modification should be carried out by rotating the trial body and/or replacing adaptors and trial heads until the best arrangement is obtained.

Before removing all the trial components, if the humeral body position has been modified, mark the new position using the electric scalpel near one of the existing marks and memorize the body insertion depth into the humeral canal.

Remove all the trial components.

#### **5. INSERTION OF THE FINNED STEM**

From the sterile box, take out the final finned stem of the size corresponding

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[16]

[17]

to the last trial stem used.  
Open also the package of the trauma humeral body of the required height.

### 5.1. Insertion of the final stem and body

Assemble the final humeral body on the finned stem. By forcibly tightening the impactor, the Morse taper is firmly fixed. Remove the impactor and tighten the safety screw.

Assemble the complete implant onto the prosthesis introducer (A2) screwing the proximal screw which locks the system (Fig. 17).

Before inserting the implant, assemble the alignment rod on the introducer handle. Place the alignment rod parallel to the forearm flexed at 90°. In this position the implant will have a 30° retroversion.

If an inferior retroversion is required, the rod must be extra-rotated (Fig. 18).

Impact the stem into the canal. Insert the gauge in the corresponding seat of the prosthesis introducer. Continue the implant insertion until the gauge reaches the desired height (Fig. 19).

Make reference to the marked point on the humeral stump to determine the insertion height.

Unscrew the proximal screw and remove the prosthesis introducer.

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[19]

[20]

## 6. INSERTION OF THE ADAPTOR TAPER IN THE HUMERAL HEAD

From the sterile box, take out the final adaptor taper and head of the type and size selected during the trial reduction.

Insert the adaptor taper into the head (Figures 20-21).

If an eccentric adaptor taper is used, align the mark present on the adaptor taper with that of the humeral head. The neutral adaptor taper is not provided with the above mark.

Place the system in the **humeral head press (S2)** and forcibly rotate the press lever. In this way the head-adaptor taper assembly is permanently fixed (Fig. 22).

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[23]

## 7. INSERTION OF THE HUMERAL HEAD

Insert the selected humeral head in the cephalic cone (Fig. 23) placing the adaptor taper peg in the same position of the trial peg.

Make sure that the contact surfaces are well cleaned and that the head or the adaptor is not jamming against the bone, otherwise the stability of the Morse coupling may be hindered. Finally, fix the tapered coupling by gently impacting with the **humeral head beater (R2)** (Fig. 24).

**Warning!** If the cemented stem is used, insert and impact the head on the taper of the humeral body only after the acrylic cement has completely polymerised.

[24]

## 8. REDUCTION AND SUTURE

The humeral component is reduced by a gentle internal rotation and traction movement.

Reconstruction of the planes. The suture of the capsule is not carried out to avoid rigidity and limitation of the external rotation. The subscapularis muscle is carefully reinserted. The muscle is fixed to the humeral metaphysis by nonabsorbable transosseous sutures.

The use of suction drainages is recommended.

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### 9. USING THE HALF-KNURLED STEM

The insertion of the half-knurled stem (Fig. 25) requires the use of a **wrench for Knurled Stems (B2)**. The above wrench has an hexagonal seat which houses the stem proximal portion.

The stem is inserted with a rotary motion. Next to the point where the knurling is no more visible inside the canal, insert the gauge (which will indicate the correct stem insertion depth into the canal) on the wrench.

Remove the stem wrench.

It is possible to check the insertion depth inserting the trial body into the final stem: insert the trial head and adaptor in the trial body and carry out the trial reduction.

After having established the correct height of the trauma humeral body and its retroversion, open the packaging the final components.

From the sterile box, take out the final humeral body (short, medium or long).

Impact it on the half-knurled stem already inserted, paying attention to the correct retroversion. Fix the two elements with the safety screw, using the body stopper (I2) inserted in the cephalic cone of the final body.

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[1]

### 1. PRE-SURGERY PLANNING

In order to obtain the best results it is highly recommended to make a pre-surgery planning by using appropriate x-ray templates. The templates give a 5% enlargement of profiles (Fig. 1). Make good quality frontal x-rays with a proper contrast. The size of the x-rays must contain the entire length of the pre-surgery template stems. Establish the stem size and the humeral head resection level. The resection level will be used as a reference for the height of the final implant.

### 2. SURGICAL FIELD PREPARATION

#### 2.1. Positioning

The patient is placed supine in the semi-sitting (beach chair) position (Fig. 2). The patient's arm must protrude out of the surgical table to allow a good mobility in all the planes. The patient's head must be fixed rigidly.

#### 2.2. Access routes

Anterior, deltopectoral approach. The incision must be rather long, from the edge of the clavicle above the caracoid process, diagonally down to the humeral insertion of the deltoid (Fig. 2). The deltopectoral groove is localised and the cephalic vein isolated (moved

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[3]

medially or ligated and then dissected). After having incised the calvipectoral fascia, the conjoined tendon of the short head of the biceps and the coracobrachialis is moved medially.

### 2.3. Humeral head dislocation

A vertical incision of the subscapularis muscle is made one centimetre from its insertion on the small tuberosity and the muscle is tilted medially. A vertical anterior capsulotomy is performed (attention must be paid to the axillary nerve). The humeral head is frontally dislocated by external rotation and extension proceeding with caution (the humerus bone can be osteoporotic).

## 3. PERFORATION OF THE HUMERAL CANAL AND RESECTION OF THE HEAD

### 3.1. Perforation

Perforate the humerus proximal end with a pointed instrument or an osteotome. Then insert the **reamer (A1)** together with the **manual snap wrench (D1)** (Fig. 3).

Insert the reamer in the canal with a rotary motion until the correct fit is obtained (the reamer cutting edges are no more visible).

Insert the **right or left resection mask (O2)** on the reamer shaft. Then insert the **alignment rod (P2)** on the mask to determine the retroversion (Fig. 4).

**Option 1:** in case of thin humeral canals (proximal portion <16 mm),

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the reamer cannot be completely inserted. Remove the reamer and screw the **conical trial stem dia. 14 mm (B1)** on the **stem impactor (C1)**. Impact the stem into the canal, checking the insertion depth with the right or left resection mask mounted on the impactor rod (Fig. 4).

**Option 2:** in case of large dimension humerus, the reamer may take an incorrect angle in the canal. This can be avoided by screwing a trial stem having a diameter larger than that of the reamer (16 mm) on the stem impactor (C1) and impacting it in the canal. Check the obtained insertion depth inserting the **resection mask (O2)** in the impactor shaft (Fig. 4). In case of excessive insertion depth, repeat the operation with a larger size stem.

### 3.2. Retroversion

Maintaining the forearm flexed at 90°, rotate the resection mask until the forearm and the alignment rod are parallel. In this position, a resection with 30° retroversion will be carried out. If a smaller retroversion is required, the rod must be extrarotated.

### 3.3. Resection of the head

Insert the reamer with a rotary motion until the oscillating blade resting plane is approximately at the anatomic neck level. Fix the mask with the Kirschner wires (dia. 2.5 mm), not included in the instrument set (Fig. 5).

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[5]

[6]

Proceed with the humeral head resection using a thin blade placed on the upper surface of the resection mask. Perform a partial osteotomy.


Remove the Kirschner wires, the resection mask and the reamer. Then proceed with the complete resection of the head (Fig. 6).

**Warning!** Having chosen the resection with cutting mask assembled on the stem impactor, due to the partial osteotomy, remove only the impactor leaving the trial stem in place.

If an arthroprosthesis is being performed, proceed with the glenoid implant.

Make reference to the chapter "SMR GLENOID".

The SMR instrument set includes a **cover (B3)** to be inserted in the humerus in order to coat the resected surface.

[7]

#### **4. HUMERUS PREPARATION**

##### **Stem size selection**

If the glenoid has been implanted, remove the glenoid retractor and expose again the humerus. Remove the protection cover and proceed with the humeral preparation.

Screw the stem impactor (C1) on the conical trial stem (B1) (Fig. 7), impact the assembly into the canal and make sure that the **45° stop guide (N2)**, placed in the corresponding seat ob-

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
tained in the impactor shaft, is in contact with the resected surface (Fig. 8).

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Remove the stop guide and impact again the stem. If the latter sinks, insert progressively the larger stem diameters until a stem which does not sink beyond the established level is found.

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Remove the stop guide and unscrew the stem impactor, leaving the last trial stem used in place (Fig. 9).

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The previously selected retroversion is obtained using the **body aligner (I2)** placed on the resected humeral surface, by inserting the alignment rod in the corresponding hole (right or left humerus).

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If a retroversion different from 30° is required, do not insert the rod, but maintain the aligner rod adherent to the resected humeral surface.

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Then tighten the screw of the trial humeral body with the **allen wrench (M2)** (Fig. 11).

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**Option.** It is possible to assemble the stem to the body outside the surgical field, as described in chapter "SMR TRAUMA – Insertion of the trial implant (option 2)".

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
Drive the **trial adaptors, STD** (neutral

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**[12]**

or eccentric) or **long (E2)** (neutral or eccentric) type, into the **trial head (C2)** by hand (Fig. 12) and reduce the implant checking the contact with the glenoid and the articular tension (Fig. 13).

Use the **pliers for trial adaptor (F2)** to remove the trial adaptors (Fig. 14).

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### **5.2. Trial reductions and modification of the components**

Any modification should be carried out by replacing the adaptors and/or trial heads until the best arrangement is obtained (congruency of the articulation).

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**[13]**

Decide whether it is necessary to increase/decrease the retroversion by rotating the trial body. Any variation carried out should be marked on the resection plane using the electric scalpel. This reference mark will be used during the insertion of the final components.

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**Note.** The change in retroversion of the humeral component implies the redefinition of the osteotomy plane.

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**[14]**

If an eccentric adaptor is used, mark the head position on the bone near the arrow present on the head, using the electric scalpel.

This procedure allows the final head to be placed in the same position.

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[16]

## 6. INSERTION OF THE FINNED STEM

From the sterile box, take out the final finned stem of the size corresponding to the last trial stem used. From the box, take out also the final finned humeral body.

### Assembling the final components

Outside the surgical field, assemble the humeral body and the final finned stem (Fig. 15). By forcibly tightening the stem impactor, the stem and the humeral body are permanently assembled (Fig. 16). Remove the impactor and tighten the safety screw using the allen wrench (Fig. 17).

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[18]


[19]

Assemble the **prosthesis introducer (A2)** (Fig. 18) onto the prosthesis by screwing the internal proximal screw with the hexagonal wrench (Fig. 19).

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Impact the implant into the humeral canal. To carry out a further assessment of the retroversion and ensure that the implant has been correctly inserted, apply the alignment rod to the handle of the prosthesis introducer, placing the rod parallel to the forearm (Fig. 20). If an inferior retroversion is required, the rod should be extra-rotated. Continue the insertion until the prosthesis introducer stop is in contact with the resected humeral surface (Fig. 21).

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[21]

[22]

## 7. INSERTION OF THE ADAPTOR TAPER IN THE HUMERAL HEAD

From the sterile box, take out the final adaptor taper and head of the type and size selected during the trial reduction.

Insert the adaptor taper into the head (Figures 22-23).

If an eccentric adaptor taper is used, align the mark present on the adaptor taper with that of the humeral head. The neutral adaptor taper is not provided with the above mark.

Place the system in the **humeral head press (S2)** and forcibly rotate the clamp lever. In this way the head-adaptor taper assembly is permanently fixed (Fig. 24).

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[25]

## 8. INSERTION OF THE HUMERAL HEAD

Insert the selected humeral head on the cephalic cone (Fig. 25). If an eccentric taper is used, place the head so that the chamfer of the prosthetic humeral head is placed on the previously marked reference.

Make sure that the contact surfaces are well cleaned and that the head or the adaptor is not jamming against the bone, otherwise the stability of the Morse coupling may be hindered.

Finally, fix the tapered coupling by gently impacting with the **humeral head beater (R2)** in the direction of the head axis (Fig. 26). Then impact in the direction of the diaphyseal axis to fill the small gap present between the head and the humeral resection.

**Note.** If the stem does not sink further, impaction should not continue.

## 9. REDUCTION AND SUTURE

The humeral component is reduced by a gentle internal rotation and traction movement.

Reconstruction of the planes. The suture of the capsule is not carried out to avoid rigidity and limitation of the external rotation. The subscapularis muscle is carefully reinserted. The muscle is fixed to the humeral metaphysis by nonabsorbable transosseous sutures. The use of suction drainages is recommended.

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## **10. USING THE CEMENTED STEM**

If cemented stems are used (Fig. 27), the surgical technique presents few and simple variations. Proceed as described above. In section 6, replace “finned stem” with “cemented stem”. Use a stem with smaller diameter than the trial one.

Instruments and methods are the same.

The canal is filled with cement, then the assembled components are inserted with the prosthesis introducer. Proceed with the surgical technique only when the cement has completely polymerised.

Attention must be paid to the humeral body proximal cleanliness. The seat for the adaptor taper must remain without cement.

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## “SMR Glenoid” Surgical Technique

## Tecnica operatoria “SMR glenoide”

The glenoid time of a shoulder total prosthetic implant is performed after having resected the humeral head and, in any case, always before having inserted the final humeral prosthesis. As for the access route and the preparatory stages, refer to the SMR Elective surgical technique.

However, consider that the interposition of the glenoid prosthesis between the scapula and the humeral head necessarily implies a reduction of the articular space. Therefore, a larger resection of the humeral head may be necessary to avoid an excessive tension of the glenohumeral articulation. The humeral head larger resection requires the utmost attention in order to preserve the infraspinatus posterior insertion.

### 1. EXPOSURE

Expose the glenoid with the appropriate retractors.

Using the **Fukuda retractor (F3)** it is possible to maintain the proximal portion of the humerus outside the articular space, thus maintaining the best exposure of the glenoid surface.

The retractor is generally fixed to the posterior rim of the glenoid neck in case of anterior surgical approach (deltopectoral), or to its lower rim, in case of superolateral approach (transdeltoid).

The SMR instrument set includes a **humeral cover (B3)** which coats the resected humeral surface. This

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cover allows to lever the humerus using the Fukuda retractor, without damaging it. It is advisable to remove the glenoid labrum to obtain a good exposure of the articular surface.

## **2. CEMENTLESS GLENOID IMPLANT**

The SMR cementless implant is recommended in cases of good bone trophism, sufficient to ensure an adequate tightening of the additional screws.

### **2.1. Determination of the glenoid centre**

After having completely exposed the glenoid, mark two cross orthogonal lines along the main axes (longitudinal and transversal) using the electric scalpel (Fig. 1).

This operation does not require special measurements as its purpose is to determine approximately the glenoid centre and thus avoid errors during the implantation.

[1]

**Note.** Make sure to remove any osteophytes present on the bone contour.

### **2.2. Insertion of the guide wire**

To prepare the glenoid seat, use a 15 cm long guide wire having a diameter of 2.5 mm, not included in the instrument set.

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[2]

Place the wire in a slightly decentralised position (Fig. 2) and insert it orthogonally to the glenoid surface for about 3 cm (Fig. 3).

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
[3]

The direction taken by the guide wire determines the subsequent version taken by the glenoid component (Metal Back). Therefore, it is advisable to perform a previous computerised axial tomography (CAT) to assess any deformations of the articular surface due to degenerative pathologies or trauma.

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Any corrections should be carried out during this surgical time as the implant technique does not allow any subsequent correction.

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[4]

Option: for glenoids having an upward vertical slope, correct the anatomic variant by slightly bending the wire caudally before inserting it in the bone (Fig. 4).

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[5]

After having fixed the Kirschner wire, insert the **Small or STD glenoid reamer (G3)** in the wire, depending on the glenoid dimensions (Fig. 5).

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
[6]

Apply the **glenoid reamer shaft (H3)** to the reamer and proceed by reaming the bony surface (Fig. 6). This operation is meant to remove the cartilage and expose the subchondral bone.

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[7]

Remove the reamer used and then carry out the core using the **Small-R glenoid drill (D3)**, inserted in the previously implanted guide wire (Fig. 7). Proceed with the coring operation until reaching the stop position. If the Small-R (M-B) Metal Back component is implanted, remove the guide wire and proceed with the sequence. If a component with a larger size peg (Small size M-B) or a component with larger peg and surface (STD size M-B) is going to be implanted, insert the **glenoid drill (D3)** (Fig. 8) in the wire and widen the hole.

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[8]

Glenoid Drill  
Carotatore glenoide  
cod. 9013.15.120

#### **2.4. Insertion of the cementless glenoid**

From the sterile box, take out the M-B component corresponding to the required size and insert it on the correspondent **M-B glenoid impactor (C3)** (Fig. 9).

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
[9]

**Warning.** The glenoid prosthesis cannot be hooked to the impactor. Therefore, pay the utmost attention during this surgical time.

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[10]

Fit the central peg of the prosthesis inside the hole previously made, gently impacting with a hammer on the positioner handle (Fig. 10).

The longer axis of the prosthesis must coincide with the longer axis of the glenoid.

**Note.** The Metal Back glenoid prosthesis is covered with hydroxyapatite. Therefore, reduce to a minimum the direct contact with hands or metal surfaces.

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**2.5. Insertion of the fixing screws**

After having positioned the cementless glenoid, drill the seats for the fixing screws using the flexible mandrel (K3) with the 3.5 diameter helix drill (M3) inserted in the SMR drill guide (I3) (Fig. 11).

The spherical seat for the heads of the screws on the metal shell allows the implant direction to be chosen within a certain angular range. However, it is advisable to maintain a direction approximately orthogonal to the metal shell.

After having prepared the seat for the first screw, insert the screw with the screwdriver (L3) without tightening it completely (Fig. 12). Then prepare the second seat and insert the screw.

The tightening of the screws must be carried out at the same time to ensure the best support of the metal glenoid on the previously reamed bone seat.

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
[12]

### **2.6. Insertion of the UHMWPE liner**

From the sterile box, take out the articular liner having the same size of the already implanted cementless glenoid and, after having carefully cleaned the internal part and the shell rim from fat and soft tissues, insert the liner pressing it with the thumb until the noise of a snap is heard (Figures 13-14).

**Note.** The liner insertion procedure is not reversible. Should it be necessary to remove intraoperatively the liner, open the package of a UHMWPE liner of the same size and carry out a new implantation.

During the following humerus stages, attention should be paid not to damage the liner with metal objects. Then remove the retractors. The Fukuda retractor must be removed carefully to avoid any sticking with the glenoid prosthesis.

### **3. CEMENTED GLENOID IMPLANT**

The cemented implant is recommended in cases of scarce bone trophism, i.e. when an adequate tightening of the central peg and the additional screws of the SMR cementless model cannot be ensured.

The procedure for preparing the glenoid presents few and simple variations. Proceed as described for the cementless surgical technique up to section "2.3 Glenoid surface preparation".

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### Choosing and cementing the prosthesis

From the sterile box, take out the UHMWPE glenoid of the size corresponding to the reamer used. Make sure that the prosthesis is in complete contact with the prepared seat. Place the acrylic cement in the prepared cavity and then insert the final glenoid component maintaining it pressed with the **cemented glenoid pusher (E3)** until the cement has completely polymerised.

Then remove the retractors. The Fukuda refractor must be removed carefully to avoid any sticking with the glenoid prosthesis.

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
[1]

### 1. PRE-SURGERY PLANNING

In order to obtain the best results it is highly recommended to make a pre-surgery planning by using appropriate x-ray templates. The templates give a 5% enlargement of profiles (Fig. 1).

Make good quality frontal x-rays with a proper contrast. The size of the x-rays must contain the entire length of the pre-surgery template stems.

Establish the stem size and the humeral head resection level. The resection level will be used as a reference for the height of the final implant.

### 2. SURGICAL FIELD PREPARATION

#### 2.1. Positioning

The patient is placed supine in the semi-sitting (beach chair) position (Fig. 2). The patient's arm must protrude out of the surgical table to allow a good mobility in all the planes. The patient's head must be fixed rigidly.

#### 2.2. Access routes

The SMR reverse prosthesis can be implanted following any of the usual approaches:

- deltopectoral
- transdeltoid.

Each one presents advantages and disadvantages and the choice is left to the surgeon. The surgical technique described below is independent from the chosen approach.

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
### 3. PERFORATION OF THE HUMERAL CANAL AND RESECTION OF THE HEAD

#### 3.1. Perforation

Perforate the humerus proximal end with a pointed instrument or a small osteotome. Then insert the **reamer (A1)** together with the **manual snap wrench (D1)** (Fig. 3).

**Note.** In case of scarce trophism, the perforation can be made directly with the reamer.

Insert the reamer in the canal with a rotary motion until the correct fit is obtained (the cutting edges are no more visible).

Insert the **resection mask for reverse prosth (H4)** on the reamer shaft. Then insert the **alignment rod (F4)** in the latter to determine the retroversion (Fig. 4).

**Option 1:** in case of thin humeral canals (proximal portion <16 mm), the reamer cannot be completely inserted. Remove the reamer and screw the **14 mm conical trial stem (B1)** on the **stem impactor (C1)**. Impact the stem into the canal, checking the insertion depth using the resection mask (H4) inserted in the impactor shaft (Fig. 4).

**Option 2:** in case of large dimension humerus, the reamer may take an incorrect angle in the canal. This can be avoided by screwing a trial stem having a diameter larger than that

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[3]

of the reamer (16 mm) on the stem impactor (C1) and impacting it in the canal. Check the obtained insertion depth inserting the resection mask (H4) in the impactor shaft (Fig. 4).

In case of excessive insertion depth, repeat the operation with a larger size stem.

### 3.2. Retroversion

Insert the alignment rod in the hole marked with 0° or 20°. With the forearm flexed at 90°, rotate the resection mask until the horizontal rod is parallel to the forearm.

In this position the retroversion of the humeral resection has a 0° or 20° angle, depending on the selected hole.

**Note.** The 0° or 20° retroversion angle is indicative and the surgeon can change it as required.

### 3.3. Head resection

Insert the reamer with a rotary motion until the oscillating blade resting plane is approximately at the anatomic neck level. Fix the mask with Kirschner wires (dia. 2 mm), not included in the instrument set (Fig. 4).

Proceed with the humeral head resection using a thin blade placed on the upper surface of the resection mask. Perform a partial osteotomy. Remove the Kirschner wires, the resection mask and the reamer. Then proceed with the complete resection of the head (Fig. 5).

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**Warning!** Having chosen the resection with cutting mask assembled on the stem impactor, due to the partial osteotomy, remove only the impactor leaving the trial stem in place.

**Note.** In case of doubts, especially in presence of tight articulations, even after an adequate release, it is advisable to carry out a large resection.

The glenoid is then implanted. The SMR instrument set includes a **humeral cover (B3)** to be inserted in the humerus in order to coat the resected surface.

#### **4. GLENOID IMPLANT**

##### **4.1. Insertion of the cementless glenoid**

In the reverse prosthesis, the glenoid component uses the cementless implant procedure already indicated in the chapter "SMR GLENOID". Therefore, follow the procedure for the Metal Back component described in the chapter "SMR GLENOID" up to the section "2.5 - Insertion of the fixing screws".

##### **4.2. Removal of the bone contour**

Hook the trial glenosphere by screwing the **glensphere impactor-extractor (T4)** on one of the holes made on the outer rim and position it on the previously implanted Metal Back glenoid (Fig. 6).

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
[6]

[7]

Insert the **guide screw for trial glenosphere (R4)**, (P4 for Small-R Metal Back) through the polar hole and tighten completely using the screwdriver (Fig. 7). In this configuration the glenosphere is free to slide along the screw and can be used to ream the glenoid surface.

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## 5. HUMERUS PREPARATION

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### 5.1. Stem size selection

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After having removed the glenoid retractor, the humerus is again exposed. Remove the protection cover and proceed with the humeral preparation. Screw the stem impactor (C1) on the conical trial stem (B1) (Fig. 9), impact the assembly with a beater and make sure that the **30° stop guide (N4)** placed in the corresponding seat obtained in the impactor shaft is in contact with the resected surface (Fig. 10).

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[9]

[10]

Remove the stop guide and impact again the stem. If it sinks, use the larger trial stem diameters until a

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[11]

stem which does not sink beyond the established level is found.

**5.2. Preparation of the seat for the insertion of the reverse body**

Unscrew the stem impactor, leaving the last trial stem used in place.

From the Reverse instrument set, take out the lower part of the **humeral body impactor (L4)** and remove its movable plastic part with a translation movement of the anterior rim and traction along the axis (Fig. 11).

Place the beater on the depression obtained on the **guide for conical reamer (I4)**, and press the proximal end with the forefinger until a click is heard (indicating the pin engagement caused by the spring grains) (Fig. 12).

[12]

Place the impactor near the stem housed in the canal and fix the modules beating with an hammer from above (Fig. 13). Use the external rod to determine the retroversion.

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

The positioner is removed by slightly bending the handle laterally (thus eliminating the constraint reaction due to the spring grain) and lifting it (Fig. 14).

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Mill using the **conical reamer (J4)** mounted on the manual snap wrench (Fig. 15).

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It is recommended to carry out this operation gently by hand, without using motorised manipulators.

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**[14]**

Make sure that the humeral seat is correctly prepared by inserting a Kirschner wire (2.5 mm max. diameter) in the reamer clearance hole, near the cutting edges (Fig. 16). The passage of the Kirschner wire means that the humeral seat has not been completed.

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[17]

[18]

After completing the humeral metaphysis milling, remove the guide for conical reamer. First screw in the **extractor for humeral body (C4)** and then the **universal stem for extractor (D4)**, paying attention not to produce torsions on the humerus (Fig. 17). When the grip is perceived, force the screwing until the component disengages from the inferior stem (Fig. 18). Therefore proceed with the piece removal (Fig. 19).

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## 6. REDUCTION OF THE TRIAL IMPLANT

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### 6.1. Insertion of the final glenosphere

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Insert the peg in the seat obtained in the glenosphere and tightly fix the components screwing the knurled portion of the glenosphere impactor-extractor (T4) until reaching the end of stroke (Fig. 20).

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[20]

[21]

Insert the assembled system inside the glenoid cavity and impact with a beater to favour the Morse taper coupling (Fig. 21).

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[22]

Unscrew the impactor and check by hand the stability of the Morse taper. To increase the solidity of the system, insert the final safety screw contained in the box (Fig. 22).

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### 6.2. Insertion of the trial reverse humeral body and liner

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[23]

[24]

Take the trial reverse humeral body (G4) from the instrument set and, using the impactor (lower part) gently impact into the canal (Fig. 23).

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Insert the trial STD liner (O4) (Fig. 24) by hand and carry out the trial reductions on the previously implanted final glenosphere (Fig. 25).

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[25]

First assess the articular tension and correct any relaxation by replacing the STD size with a +3 or +6 size. To remove and replace the trial liner, lever with a pointed object through one of the two holes or use a chisel placed on the lateral rim. At the same time, assess the possible notch on the medial rim and, if necessary, rotate the liner until placing the chamfered side on the point of contact with the bone (ROM increase). Then assess the articular congruence and, in case of instability, rotate the liner positioning the chamfered side near the large tuberosity.

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removed, use the extractor (C4) as described for the guide for conical reamer.

**Note.** Before removing the trial body, mark the position of the main groove on the bone in order to reproduce the same position for the final one.

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## 7. USING THE ECCENTRIC GLENOSPHERE

### 7.1 Trial component implant

Hook the trial eccentric glenosphere by screwing the glenosphere impactor-extractor (T4) on the hole made on the upper rim and place it near the Metal Back glenoid.

Insert the guide screw for trial glenosphere (R4), (P4 for Small-R Metal Back), through the hole and tighten completely maintaining the eccentricity of the component in a glenoid lower quadrant (Fig. 26).

### 7.2. Trial reduction

Insert the trial reverse body and liner as indicated and carry out the reduction of the new articulation. Before carrying out this operation, make sure that the sphere has been correctly fixed on the Metal Back glenoid.

In case of instability or notch of the liner on the scapula, dislocate and find a better position for the sphere on the Metal Back.

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
[26]

[27]

found, before removing the component, fix a landmark near one of the two contour steps (Fig. 27).

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**7.3. Insertion of the final eccentric glenosphere**

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
Open the sterile box and assemble the peg to the sphere as described for the standard final glenosphere.

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Then place the system near the Metal Back glenoid making reference to the landmark previously marked and impact in order to favour the coupling (Fig. 28).

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[28]

Unscrew the impactor and check by hand the stability of the Morse taper. To increase the solidity of the system, insert the final safety screw contained in the box.

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**7.4. Removal of the trial implant**

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Remove the plastic liner and screw in the stem impactor (C1) (Fig. 29). Extract the implant by percussion.

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**8.1. Assembling the final stem and reverse humeral body**

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Outside the surgical field, assemble the final humeral body to the final finned stem using the whole body impactor (L4). By forcibly tightening the positioner, the stem and the humeral body are permanently assembled (Fig. 30).

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**[31]**

Then slightly loosen the tightening. Thanks to this simple operation it is possible to avoid dangerous torsions which may damage the humerus while the impactor is removed.

**8.2. Insertion of the final humeral components**

Impact the stem and the body into the canal. The retroversion is determined by placing the body main notch on the previously marked bone landmark or using the alignment rod.

After reaching the complete insertion depth of the reverse body on the resection plane, unscrew and remove the impactor handle, insert the reverse body safety screw and tighten using the appropriate hexagonal wrench (Fig. 31).

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**[32]**

**8.3. Insertion of the final reverse liner**

Using the fingers, insert the selected polyethylene liner in the reverse humeral body (Fig. 32). The body internal part must be clean and dry.

The liner groove must be placed on the body notch used during the reduction. This allows the reverse liner to be inserted with the same phase selected during the trial reduction.

Exert a slight pressure and make sure that the liner is correctly positioned. Then impact the liner with the **cemented glenoid pusher (E3)** (Fig. 33).

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

**[34]**

#### **8.4. Reduction**

Carry out the normal articular movements, carefully evaluating the impingement of the humerus against the acromion (Fig. 34).

In this case, carry out an acromioplastic surgery, if necessary.

It is also necessary to check the presence of humeral or glenoid osteophytes, which must be removed.

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[1]

### FOREWORD

The prosthetic treatment with CTA heads is clinically indicated for the treatment of cuff arthropathies in which it is not possible to implant a reverse prosthesis due to anatomic anomalies or scarce bone trophism of the glenoid.

### 1. FIRST IMPLANT PRE-SURGERY PLANNING

In order to obtain the best results it is highly recommended to make a pre-surgery planning by using appropriate x-ray templates. The templates give a 5% enlargement of profiles (Fig. 1).

Make good quality frontal x-rays with a proper contrast. The size of the x-rays must contain the entire length of the pre-surgery template stems.

Establish the stem size and the humeral head resection level. The resection level will be used as a reference for the height of the final implant.

### 2. ACCESS ROUTES AND EXPOSURE

Anterior, deltopectoral approach. The incision starts from the anterior edge of the clavicle above the coracoid process and continues diagonally down to the humeral insertion of the deltoid. The exposure of the humeral head is obtained with a cautious extrarotation and extension of the limb.

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
The traction and repositioning of the elbow may favour the dislocation of the cephalic component, which has generally migrated upward and is in contact with the subacromial arch.

### 3. CANAL PERFORATION AND DEFINITION OF THE OSTEOTOMY PLANE

Perforate the humerus proximal end using the reamer (A1) near the large tuberosity-head middle line (Fig. 2) and assess the insertion depth reached with the help of the resection mask (O2).

Insert the reamer so that the mask resting plane is near to the anatomic neck. Fix the mask and start the bone resection as described in chapter "SMR ELECTIVE" (Fig. 3).

Then proceed with the resection of the head as described in the section with the same title.

### 4. HUMERUS PREPARATION

#### 4.1. Insertion of the trial implant and reduction

After having completed the osteotomy, select the stem size as indicated in section "SMR ELECTIVE – Stem size selection", then assemble a medium size trial humeral body (G2) to it. Then impact it into the canal (Figures 4-5).

Using the head gauge (D2), measure the size of the cephalic component

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[6]

and, using the **pliers for trial adaptors (G2)**, insert a **neutral STD adaptor (E2)** in the corresponding **trial head (C2)**. Assemble the component in the bevel of the humeral body and carry out a trial reduction.

Assess the linear excursion of the new articulation during the limb traction and release movement and correct any articular relaxation using a long adaptor. If this is not sufficient, decide whether it is necessary to use an eccentric adaptor in order to lift or lower the head rotation centre along the craniocaudal axis.

If necessary, correct the consequent offset variation with a long or short trial body. Finally check the head measure carrying out intra-extra rotation movements in relation to the overall space filled by the component and the lengthening undergone by residual tendons.

#### **4.2. Definition of the osteotomy lateral plane**

After having determined the required size of the head and the adaptor taper, keep the components in the chosen position and carry out a raiser for the cut of the large tuberosity using the oscillating blade through the slot obtained on the trial head (Fig. 6).

Then remove the component and complete the bone resection. After having removed the trial components, the bone metaphysis appears as in figure 7.

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[8]

## 5. REDUCTION OF THE FINAL IMPLANT

### 5.1. Assembling and inserting the stem and the humeral body

From the sterile boxes, take out the stem and the Trauma humeral body of the required sizes. Assemble the parts as described in section “SMR TRAUMA – 5.1. Insertion of the final stem and body” and place the assembly into the canal following the standard procedure (Fig. 8).

Insert the component down to the level of the osteotomy anatomic plane that has been carried out.

[9]

### 5.2. Assembling the adaptor taper and the head

From the box, take out the adaptor taper and the head of the required sizes and tightly fix the components using the **humeral head press (S2)** (Fig. 9).

If an eccentric adaptor taper is used, align the mark with the lateral middle line to obtain a cranial positioning of the head, or medially to obtain a caudal positioning (Fig. 10).

craniale

caudale

[10]

### 5.3. Humeral head insertion

Insert the selected humeral head in the cephalic cone.

If a eccentric adaptor is used, place the head taking care to respect the required offset.

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[11]

Make sure that the coupling surfaces are well cleaned and then fix the coupling by gently impacting with the **humeral head beater (R2)** (Fig. 11). Fill the gap, if any, further impacting in the diaphyseal axis direction.

## 6. REVISION WITH CTA HEADS

### 6.1. Removal of the glenoid implant

Expose the glenoid with appropriate retractors and remove the mobilised components (glenosphere, Metal Back and screws). In so far as the conditions of the residual bony portion allow, carry out a bone graft before proceeding with the humeral implant.

### 6.2. Removal of the humeral liner

Place the chisel on the superior rim of the polyethylene liner and, using a hammer, impact until the component is removed from the housing of the reverse body (Fig. 12).

Carefully clean the body housing and assess the height of the implant in relation to the residual bony portion of the glenoid and the subacromial arch.

### 6.3. Insertion of the humeral head

Carry out the x-rays of the contralateral articulation and with the help of the pre-surgery templates, assess the required size of the cephalic component (head).

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[13]

From the box, take out the reverse body 36 mm adaptor, screw the knurled part of the **extractor for humeral body (K2)** onto the adaptor and impact gently (Fig. 13).

Then take the CTA head of the required size and place it on the adaptor, making sure that the present bone outer rim is not jamming with the head lateral extension. If this is the case, remove part of the bony tissue with a Luer.

Then impact the head using the humeral head beater (R2) (Fig. 14). In the revision implant with CTA head it is not possible to change the eccentricity of the cephalic component with respect to the cervico-diaphyseal axis.

**Note.** In case of excessive articular relaxation, before implanting the 36 mm adaptor, apply a +9 mm humeral extension using the method described for the 36 mm adaptor.

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“SMR Revision and Tumoral”  
Surgical Technique

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[1]

**FOREWORD**

The technique described in this section deals with the use of the instrument set 9013.60.000 for the revision of Lima SMR shoulder implant or other prosthetic implant.

Therefore the procedure described implies the previous removal of the existing implant and the implantation of a SMR modular prosthesis with finned (Fig.1) or cemented revision stem with the coupling of modular spacers in case of large bone resection due to tumour.

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[3]

## 1. PRE-SURGERY PLANNING

In order to obtain the best results it is highly recommended to make a pre-surgery planning to assess the component sizes necessary for the prosthesis implantation (Fig. 3) and, in case of reverse prosthesis implantation, the height of the humeral implant in relation to the glenoid.

## 2. REVISION WITH REVERSE PROSTHESIS

### 2.1. Cementless stem size selection

Screw the trial revision stem (D6) selected during the pre-surgery planning on the stem impactor (C1) and gently impact into the canal.


Assess the component tightness with vigorous traction movements (Fig. 4). If the component finned portion remains exposed outside the canal by more than 3 cm, decide whether to use a longer or a cemented stem.

### 2.2. Glenoid implant

Expose the glenoid cavity using the appropriate retractors and proceed with the Metal Back glenoid implantation following the procedure described in chapter "SMR REVERSE – 4. GLENOID IMPLANT".

### 2.3. Insertion of the final glenosphere

Proceed by inserting the final glenosphere in the already implanted

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[5]

prosthetic glenoid. Follow the procedure indicated in section “SMR REVERSE – Insertion of the final glenosphere”.

**2.4. Assembling the trial stem and humeral body**

Outside the surgical field, assemble the selected stem to the **trial reverse body (G4)** (Morse taper coupling) and impact into the canal.

During the insertion operations, take care to reach the insertion depth which ensures the fit of the stem and use the **alignment rod (F4)** parallel to the forearm to determine the retroversion angle (Fig. 5).

[6]

Remove the impactor and use a **trial STD liner (O4)** (Fig. 6) to carry out a trial reduction (Fig. 7).

Assess the congruence and stability of the new articulation and change the size and position of the liner as required.

**2.5. Assembling and inserting the final humeral components**

Remove the trial implant using the stem impactor. Open the boxes of the stem and the reverse body, assemble the modules by forcibly tightening the **humeral body impactor (L4)**, release slightly and then impact into the canal making sure to reproduce the same implant height.

Remove the beater, position and tighten the safety screw with the **allen wrench (E4)**. Open the package of the

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
[7]

liner of the required size and complete the humeral implant following the instructions given in section “SMR REVERSE – 8.3. Insertion of the final reverse liner”.

### **2.6. Reduction of the final implant**

Carry out the normal articulation movements, carefully evaluating the impingement of the humerus against the acromial vault (acromioplasty).

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## **3. USING THE CEMENTED STEM**

### **3.1. Insertion of the trial implant**

If the humerus bone trophism does not allow the implantation of a finned stem, the modularity of the SMR system envisages the use of the cemented revision stem alone or together with a spacer.

Assemble the **trial cemented stem (C6)** and the trial reverse body and insert the assembly into the canal. Position an STD liner and carry out a reduction on the previously implanted glenosphere (Fig. 8).

The trial stem is not equipped with antirotation fins and therefore only an indication of the articular tension is obtained during this step in relation to the glenoid implant while it is not possible to verify the dynamic stability of the prosthesis.

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[9]

### 3.2. Assembling and inserting the final cemented stem

From the box, take the stem and the spacer of the selected size, slide the spacer over the stem starting from the tip and assemble the modules by forcibly tightening the system positioned inside the **press for revision stems (A6)** (Fig. 9).

Insert the humeral body and tighten the Morse taper following the standard procedure.

Then fill the canal with cement and, using the whole body impactor (L4), insert the assembly into the canal until the lower surface of the spacer is in contact with the resected surface.

Wait for some minutes to allow the cement to set and then proceed with the surgical sequence as indicated in "SMR REVERSE" until the implant is reduced (Fig. 10).

[10]

## 4. TUMORAL (CTA HEADS)

### 4.1. Stem and spacer size selection

First carry out the humerus osteotomy and assess approximately the length of the resected bony portion by comparison with the **trial spacers (B6)** (Fig. 11). This measure is useful to determine in advance the measure of the spacer to be used.

Then assemble the spacer on a trial cemented stem (C6) having a diameter of 7 mm, assemble a **trial humeral body (G2)** on the proximal portion

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[ 12 ]

and insert the assembly by hand into the canal. In case of particularly large size humerus, use a 10 mm diameter stem.

#### **4.2. Trial adaptors and heads**

Insert a **neutral or eccentric adaptor (E2)** in the trial head and assemble it to the body (Fig. 12). Reduce and check the contact with the glenoid and the subacromial arch.

If an eccentric adaptor is used, memorize the peg slot on the humeral body.

#### **4.3. Insertion of the final implant**

From the box, take the stem and the spacer of the selected size, slide the spacer along the stem starting from the tip and assemble the modules by forcibly tightening the system positioned inside the press for revision stems (A6) (Fig. 9).

Draw near the trauma humeral body and tighten the Morse taper following the standard procedure (Fig. 13).

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**[14]**

Assemble the CTA head to the selected adaptor using the humeral head press (T2), impact onto the cephalic cone and proceed with the reduction (Fig. 14).

As an alternative to the cemented stem, it is possible to use the cementless revision stem (Fig. 15).

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**[15]**

“SMR/Randelli Reverse  
Revision” Surgical Technique

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[1]

**1. REMOVAL OF THE COMPONENTS**

**1.1. Removal of the humeral head and the adaptor taper**

Insert the **head extractor (B4)** in the space below the humeral head and lever (Fig. 1). Remove the head and the adaptor taper still assembled.

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**1.2a. Removal of the Randelli humeral body**

Using the hexagonal wrench, unscrew and remove the safety screw inside the humeral body.

Screw tight the **extractor for Randelli reverse humeral body (E5)** inside the body, maintaining the direction of the stem axis (Fig. 2). Then insert the **universal stem for extractor (A5)** inside the extractor and tighten it by hand.

When the advancement starts to be difficult, insert the manual snap wrench (U4) on the hexagonal seat and, holding the humeral body extractor, force the universal stem by tightening it until the disassembly is carried out (Fig. 3).

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[3]

**1.2b. Removal of the SMR humeral body**

Using the **hexagonal wrench (E4)**, unscrew and remove the safety screw inside the humeral body.

Screw tight the **humeral body extractor (C4)** inside the humeral body, maintaining the direction of the stem axis.

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[4]

extractor (D4) inside the extractor, position manual snap wrench (U4) and tighten it.

When the advancement starts to be difficult, insert the manual snap wrench (U4) on the hexagonal seat and, holding the humeral body extractor, force the universal stem by tightening it until the disassembly is carried out (Fig. 4).

Then proceed with the reverse prosthesis implantation following the standard procedure.

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## 2. GLENOID IMPLANT

### 2.1. No previous glenoid implant

If no glenoid has been previously implanted, proceed as indicated in chapter "SMR GLENOID".

### 2.2a. Previous SMR MB glenoid implant

Remove the UHMWPE liner inserting a small osteotome between the liner and the MB glenoid (Fig. 5). Proceed with the surgical sequence.

### 2.2b. Previous Randelli MB glenoid implant

Remove the UHMWPE liner inserting a small osteotome between the liner and the MB glenoid (Fig. 5). The Randelli MB glenoid cannot house the SMR glenosphere. Use a glenosphere with an appropriate connector which can be fixed inside the Randelli glenoid. Proceed with the surgical sequence.

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[8]

### 2.3. Removal of the bone contour (SMR MB)

Hook the trial glenosphere by screwing the **glenosphere impactor-extractor (T4)** on one of the holes made on the outer rim and place it on the previously implanted Metal Back glenoid (Fig. 6).

Then insert the **guide screw for trial glenosphere (R4)**, (P4 for Small-R Metal Back) through the polar hole and tighten completely using the screwdriver (L3) (Fig. 7). In this configuration the glenosphere is free to slide along the screw and can be used to ream the glenoid surface.

Then insert the **T driver for trial glenosphere (S4)** in the appropriate housing and, with clockwise and counterclockwise smooth movements, proceed with the ablation of the outer bony portion (Fig. 8).

Remove the trial glenosphere and ablate the bone contour, if present, near the Metal Back lower rim.

## 3. HUMERUS PREPARATION

### 3.1a. SMR stem

Follow the sequence indicated in the section "SMR REVERSE- 5.2 Preparation of the seat for the insertion of the reverse body".

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[10]

Then insert the **checking block (K4)** on the reamer guide peg and place the **alignment rod (F4)** in the required seat, from 0° to 20° (Fig. 10). Maintaining the peg constrained in the position indicated in the figure, rotate the aligner until the horizontal rod is parallel to the forearm flexed at 90°. In this configuration, the retroversion of the implant will have a 0° or 20° angle, depending on the seat selected.

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Tighten the reamer guide with the screwdriver.

Mill with the proper **reamer (J4)** mounted on the manual snap key.

It is recommended to carry out this operation gently by hand, without using motorised manipulators (Fig. 11).

**Warning.** In the reverse prosthesis revision operations, the retroversion checking block is never in close contact with the resection plane.

#### **4. REDUCTION OF THE TRIAL IMPLANT**

##### **4.1a. Insertion of the final glenosphere in the SMR MB glenoid**

Open the package of the final glenosphere. The sterile box contains 3 elements:

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[12]

Insert the peg in the corresponding seat obtained in the glenosphere and tightly fix the components by screwing the knurled portion of the impactor-extractor (T4) until reaching the end of stroke (Fig. 12).

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
[13]

Insert the assembled system inside the glenoid cavity and impact with a beater to favour the Morse taper coupling (Fig. 13).

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[14]

Unscrew the impactor and check by hand the stability of the Morse taper. To increase the solidity of the system, insert the final safety screw contained in the box (Fig. 14).

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#### **4.1b. Insertion of the final glenosphere in the Randelli MB glenoid**

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- glenosphere
- connector
- safety screw.

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cod. 1374.15.010


[15]

Insert the connector of the final glenosphere (Fig. 15) in its seat inside the Randelli MB glenoid so that the two protuberances are housed in the seat of the glenoid component. Proceed with the insertion until the snap of the fit is heard.

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Then insert the final glenosphere in the peg of the connector by impacting with the **cemented glenoid pusher (E3)**. The coupling between the two elements is achieved by Morse taper.

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[16]

Then permanently fix the two elements by tightening the safety screw with the screwdriver (Fig. 16).

**4.2a. Insertion of the trial reverse humeral body and liner (SMR stem)**

Follow the indications of chapter “SMR REVERSE”.

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[17]

**4.2b. Insertion of the trial reverse humeral body and liner (Randelli stem)**

Dislocate the humerus in the upper part and insert the trial reverse humeral body (G4) (Fig. 17). The reverse body is fixed on the stem with the appropriate screw for trial body-stem (C5) using the hexagonal wrench (E4) (Fig. 18).

Insert the STD trial liner (O4) and carry out the trial reduction (Fig. 19). Assess the articular tension. In case of articular relaxation it is possible to replace the STD liner with the +3 or +6 mm size.

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[18]

**Note.** If the humeral implant is “too high” or it is impossible to reduce the implant, remove the trial body and use a short reverse body.

Before removing the trial body, trace a landmark with the electric scalpel near the lateral notch.

**4.3 Using the reverse humeral body short**

The use of a short body implies the removal of a proximal bony portion

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larger than that removed during a first implant. Extract the conical reamer and gently remove by hand more bony tissue from the humeral metaphysis. It is also possible to use an oscillating saw and take away the metaphysis, taking care to preserve any tendon insertions and the integrated stem.

#### **4.1c. Using the eccentric glenosphere (SMR MB glenoid)**

Follow the indications of chapter "SMR REVERSE" in the section with the same title.

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**[19]**

### **5. REDUCTION OF THE FINAL IMPLANT**

#### **5.1a. Insertion of the final reverse humeral body on the SMR stem**

Follow the indications of chapter "SMR REVERSE".

#### **5.1b. Insertion of the final reverse humeral body on the Randelli stem**

From the sterile box, take out the Randelli final revision reverse humeral body. Insert the body on the already implanted stem.

The body lateral notch must be placed near the landmark previously marked on the bone.

Insert the impactor (L4) inside the reverse humeral body, maintaining the same axis of the stem.

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*Procedere delicatamente, impattando*

**[20]**

with a hammer until it is completely seated on the stem (Fig. 20).

Insert the **compaction screw for Body- Stem (D5)**. Using the hexagonal wrench, tighten completely the screw (Fig. 21).

Then unscrew completely and remove the compacting screw.

This procedure ensures the correct tightening of the stem Morse taper with the humeral body.

### **5.2. Insertion of the final reverse liner**

From the sterile box, take out the reverse liner of the height previously determined with the trial reductions.

The body internal part must be clean and dry. Using the fingers, insert the selected polyethylene liner in the reverse humeral body. The liner lateral groove must be placed on the lateral of the reverse body. In this configuration the chamfered side is medially positioned (Fig. 22).

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[22]

Exert a slight pressure to make sure that the liner is correctly fixed, then impact with the cemented glenoid pusher(E3) (Fig. 23).

## 6. IMPLANT REDUCTION

Carry out the normal articular movements, carefully evaluating the impingement of the humerus against the acromion. In this case, carry out an acromioplasty surgery, if necessary.

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**Note.** In case of implant instability, dislocate, remove the liner and reposition it with the chamfered side near the large tuberosity.

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Clinical Cases

*Casi clinici*

Fig. 1  
Pre-operative x-ray  
Gender: female  
Age: 65 years old  
Pathology: arthrosis

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Fig. 2  
Post-operative x-ray

*Fig. 2  
Radiografia post-operatoria*


Fig. 3  
Pre-operative x-ray  
Gender: female  
Age: 63 years old  
Pathology: massive fracture  
of rotator cuff

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Fig. 4  
Post-operative x-ray

*Fig. 4*  
*Radiografia post-operatoria*

Fig. 5  
Pre-operative x-ray  
Gender: male  
Age: 65 years old  
Pathology: revision for massive fracture of rotator cuff

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Fig. 6  
Post-operative x-ray

*Fig. 6  
Radiografia post-operatoria*

For further information about our products, please visit our web site [www.lima.it](http://www.lima.it)

This documentation has been written by Lima-Lto Sales & Technical Department

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