MRS

Modular Resurfacing Hip Prosthesis
Protesi di rivestimento d'anca modulare

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History and Philosophy

Storia e filosofia

The MRS system was conceived and designed to provide an efficient alternative to traditional operations of total hip arthroplasty or hemi-arthroplasty in young and active patients, in an initial phase of avascular necrosis and for osteoarticular pathologies affecting cartilage.

The aim is to offer these patients the advantage of metal-metal coupling of large diameter heads with lowered and thin cups, as well as the optimal support provided by modular stems. This system allows a considerable bone saving, which is particularly useful in young and active patients who, because of their life expectancy, are likely to undergo further surgery. In this case, the conversion of the MRS in a first implant can be easily performed by resection of the femoral neck.

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



Materials

Materiali

The stem is manufactured in Ti6Al4V titanium alloy, according to strictly standardised production criteria and improved by the high-level technological know-how of LimaLto.

The external surface features a macro-roughened surface obtained through the projection of a controlled flow of corundum particles. Pressure, speed and granulometry of the projected material make it possible to obtain an even and calibrated surface texture with biologically appropriate roughness values ($Ra\ 3\div 4\mu m$) in order to support and stimulate the osteointegration.

The acetabular cup made of CoCrMo alloy, with a double coating of porous Titanium + HA with a thickness of $175\ \mu m$ of PorousTi+ $55\ \mu m$ of HA both on the cap and on the fins.

The thin cup ($2.5\ \div\ 3.5\ mm$ in thickness) has the shape of a hemisphere lowered equatorially by $3\ mm$.

Femoral heads are made of a CoCrMo alloy (compliant with ISO 5832-12 standard). The heads feature a highly polished surface finishing ($Ra\ \leq 0.02\ \mu m$), which is necessary to allow the large diameter metal-metal articular coupling they are designed for.

The metal-metal coupling has been evaluated and approved by the C.R.I.T.T.- M.T.O.S. (Centre Regional d'Innovation et de Transfert de Technologie, section Materiaux, Dépôts et Traitements de Surface) as one of the best couplings tested by them (compulsory test for registration in France).

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Biomechanical Characteristics

Caratteristiche biomeccaniche


The modular MRS system consists of three components: cervico-cephalic straight stem, resurfacing head, and acetabular cup.

The stem, available in 5 sizes, is coupled conically to the femoral head and, thanks to its macro-roughened surface, ensures an excellent biological fixation. Therefore, it is possible to extend the contact area between head and neck and increase the implant support when an ischemic event occurs in the femoral head due to the damage of extraosseous vessels during surgical treatment.

A cup with a limited thickness allows excellent saving of bone stock and optimal primary stability thanks to the interference with the acetabular seat. The double coating of porous Titanium and hydroxyapatite also ensures a proper secondary osteointegration.

Moreover, the large diameter metal-metal coupling allows an excellent proprioception by the patients and the reduction of the risk of dislocations from a surgical point of view.

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Applications

The use of a resurfacing prosthesis serves three purposes:

- ▶ maximum bone stock preservation
- ▶ anatomical dimensioning of load
- ▶ low wear

For this reason, MRS is ideal for patients suffering from an initial phase of osteoarthritis.

Contraindications include serious loco regional bone dysmorphisms, serious osteoporosis and high AVN.

Moreover, the use of a resurfacing prosthesis is not recommended in case of evident dysmetria and varus deviation of the femur neck.

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
Sizes

The stems are available in five lengths: 60, 70, 80, 90, 100mm

This way, an optimal support of the implant can be always achieved.

The heads are available in the following sizes: 42, 46, 50, 54, 58mm, each one coupling to two cup diameters.

The acetabular cups are supplied in sizes from 46 to 64mm with a 2mm pitch.

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
For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

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X Lima Deutschland GmbH