U.P.tion[®]

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
■ Pre-operative planning

Preoperative planning can easily be performed with UPTION® templates for femoral sizing. Our software BFITS can also be utilized.

It is recommended to use a radiographic marker to assess X-ray magnification on an individual basis so that the proper templates can be selected. Determine the center of rotation of the hip, and choose the UPTION® femoral stem the most suitable on the frontal and the sagittal view.

The recommended surgical approach is a posterolateral approach as it allows excellent exposure of the femoral shaft for canal preparation.

Generally, any commonly employed surgical approach can be utilized with the UPTION® stem.

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■ Removal of femoral component

Remove the femoral component in place with the standard surgical technique used by the surgeon.

Once the failed implant has been removed, the femur is cleared of any remaining cement and debris if necessary.

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■ Distal canal reaming

Flexible reamers over a guide wire are used to ream distally.

The minimum diameter of reaming is the distal diameter of the femoral stem.

The distal reaming will be made in the minimal diameter if we wish a precise preparation of the diaphyseal part.

On the contrary, if we wish to avoid any distal contact, we shall choose the diameter of maximal reaming.

The minimal and maximal diameters for each stem size are indicated in figure 1:


Figure 1

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


■ Broaching the proximal femur

The smallest broach is attached to the broach handle and the proximal femur is prepared increasing broach sizes progressively.

The final broach needs to be completely stable in the femur to ensure stem stability in the metaphysis.
Then remove broach for the trial step.

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➔ If necessary, the broach can be removed using the sliding hammer used for revision implants.
In this case, use the adaptor to connect it to the broach thread.

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■ Trial reduction

The UPTION® system provides trial implants that match the final implant dimensions.

After extracting the final broach from the femur, select the corresponding size of trial prosthesis and insert it with the broach handle.

With the proper trial stem in place, the calcar is planed flush by using the calcar reamer.


With the trial stem in place, a trial neck corresponding to the group of sizes (1-2-3, 4-5-6, 7-8-9) and a trial head are selected to determine the appropriate neck length and lateral offset.

⊕ The trial necks are magnetized on the trial stems. Each trial neck represents a group of 3 sizes with different neck lengths :

Size	Neck length head center 0	Anodisation color
1-2-3	37	Yellow
4-5-6	40.5	Green
7-8-9	44	Purple

A locating pin prevents the use of a neck which is not adapted to the trial stem in place.

⚠ **The trial neck and the calcar driver can be used only with the trial stem, not on the broach.**

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Remove the head and the trial neck and extract the trial stem with the broach handle or the slide hammer.

With the sliding hammer, use the adaptor as you would the broach.

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Clean the trial stem slot with the cleaning blade.

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

■ Femoral component insertion

Prepare and place possible grafts and particularly the calcar graft which will be consolidated under the final implant collar.

The femoral component, corresponding to the trial stem, is inserted into the femoral canal, using the stem impactor.

A trial head is placed on the neck of the stem and the hip is reduced again to assess joint stability and range of motion.

Once the trial head is removed, carefully clean the joint before placing the final femoral head on the cleaned and dried taper of the stem. Impact using the femoral head impactor.

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