
ART[®]

**Implant fémoral
anatomique**
*Anatomical femoral
implant*

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
*The proximal 1/3 of the femur is characterized by a posterior concave **sagittal curvature** in its diaphyseal part and spiral torsion around a vertical axis in its metaphyseal part. These 2 curvatures are balanced ; spiral torsion induces anteversion of the neck of the femur.*


An anatomical femoral implant must respect this architecture in order to restore complex joint mechanics and to allow physiological transmission of stresses.

*The **medial curvature** of the implant optimizes cortical contact and restores a harmonious cervico-obturator curve, thereby avoiding excessive medialization or lateralization phenomena.*

*The **cervico-diaphyseal angle of 138°** is close to the resulting force of the mechanical axis of the femur and limits varus stresses.*

*The diaphyseal part, with a tapered extremity to avoid cortical contact, presents a posterior concave **radius of curvature of 1500 mm**, reproducing the femoral diaphyseal curvature.*


*The metaphyseal part presents a **spiral torsion of 8°**; the various sections present external rotation in relation to the underlying section, resulting in a harmonious movement, reflecting anatomical reality.*

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
*The **elliptic section** of the metaphyseal of the implant gradually evolves towards a circular geometry at the tip of the stem, thereby improving fill and metaphyseal rotational stability and stress-free diaphyseal guidance.*

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Preservation of anatomy restores physiological mechanics and function, ensuring stability and security.

The ART® implant is available in **2 versions : cemented** (M30NW stainless steel) and uncemented (TA6V titanium alloy) **with total hydroxyapatite coating**. The implant can receive ceramic or cobalt-chromium Biomet France heads with a diameter 22, 28 or 32 mm.

9 sizes, right of left, for each series (cemented and full HA) of implant are available to correspond as closely as possible to anatomical variants.

To ensure a precise procedure, the instrumentation comprises trial prosthesis rasps with dimensions strictly identical to those of the implants, allowing **size-for-size insertion** of the HAP or cemented version. The automatic locking system ensures precision, security and rapidity during preparation and insertion of the implant.

Implant

La gamme		Products range	
Tige ART® à cimenter		Femoral cemented stem ART®	
Taille / Size	Gauche / Left	Droit / Right	
00	P01150P2	P01140P2	
0	P01150P1	P01140P1	
1	P0115001	P0114001	
2	P0115002	P0114002	
3	P0115003	P0114003	
4	P0115004	P0114004	
5	P0115005	P0114005	
6	P0115006	P0114006	
7	P0115007	P0114007	

La gamme		Products range	
Tige ART® HA Totale		Femoral stem ART® Full HA	
Taille / Size	Gauche / Left	Droit / Right	
00	P0119YP2	P0118YP2	
0	P0119YP1	P0118YP1	
1	P0119Y01	P0118Y01	
2	P0119Y02	P0118Y02	
3	P0119Y03	P0118Y03	
4	P0119Y04	P0118Y04	
5	P0119Y05	P0118Y05	
6	P0119Y06	P0118Y06	
7	P0119Y07	P0118Y07	

Têtes fémorales	Femoral heads	CrCo	Alumina
22,2 col -2mm	22.2 neck -2mm	P0206C22	
22,2 col 0mm	22.2 neck 0mm	P0206M22	
22,2 col +2mm	22.2 neck +2mm	P0206L22	
28 col -3,5mm	28 neck -3,5mm	P0206C28	P0209C28
28 col 0mm	28 neck 0mm	P0206M28	P0209M28
28 col +3,5mm	28 neck +3,5mm	P0206L28	P0209L28
28 col +7mm	28 neck +7mm	P0206E28	
32 col -4mm	32 neck -4mm	P0206C32	P0209C32
32 col 0mm	32 neck 0mm	P0206M32	P0209M32
32 col +4mm	32 neck +4mm	P0206L32	P0209L32
32 col +8mm	32 neck +8mm	P0206E32	

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